



Product Code	32881
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No
RDM (NSIS)	2388665
CND	Z1203020599
EAN/UPC	4719003342524

Description

MICROLIFE AFIB ADVANCED EASY (Home use)

- detection of atrial fibrillation (specific heart arrhythmia) with 1 measurement
- clinically tested according to the BHS - ISO - EHS protocols and for "special patients" such as Diabetics, Elderly, Obese, Hypothesis, Pregnant.
- MAM technology (automatic analysis of 3 consecutive measurements)
- 1 Data memory
- traffic light indicator
- low battery indicator
- mains adapter port
- world's latest technologies
- Clinically reliable measurements with AFIB & MAM
- special deflation valve allows beat detection up to 30 b/min (for bradycardic patients)

Supplied with carrying bag, preformed conical cuff 22-42 cm, 4 AAA batteries 1.5V and user manual (GB, FR, IT, DE).

Measurement range: 20-280 mmHg (pressure), 30-200 beat/minute (pulse)

Size: 160x80x32 mm

Weight: 298 g including batteries

Technical Specifications

MICROLIFE AFIB technology.

An innovative and patented technology for the early detection of Atrial fibrillation (AF).

More than 30,000 fatal strokes can be prevented within the European Union by using Microlife AFIB Technology in primary care.

The new Microlife AFIB technology allows simultaneous blood pressure measurement and AF screening.

Microlife AFIB technology detects AF with high sensitivity (97-100%) and specificity (89%) as demonstrated in several clinical studies.

Nice guidelines for "Diagnosis and assessment of hypertension" recommended to use AFIB screening for routine blood pressure measurement in primary care.

AF screening by using a Microlife AFIB device is highly recommended in patients that have to monitor their blood pressure while at home.