



Product Code	32882
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2526921
CND	Z1203020501
EAN/UPC	4719003342234

Description

MICROLIFE AFIB B3 B.P. MONITOR (Home use)

- AFIBsens technology for high accuracy in Atrial Fibrillation detection, to prevent stroke risk
- MAM (Microlife Average Mode) technology: analysis of 3 consecutive blood pressure measurements to provide accurate average result
- IHB technology for early detection of irregular heartbeat
- Gentle+ technology improves the inflation and deflation rate. Adaptive pressure control technology allows optimal speed and pressure control, ensuring comfort
- clinically validated, also for children (over 12 years), diabetics, dialysis patients and patients with end-stage renal disease (ESRD), pregnant women, elderly and obese
- conical soft cuff (22-42 cm) comfortable cuff with washable cuff cover
- 99 data memory and average display of all stored readings
- traffic light indicator for blood pressure classification

Supplied with soft conical cuff (22-42 cm), carrying bag, 4 AA batteries 1.5V and user manual (GB, FR, IT, ES, PT, DE, NL, GR, TR, FA, Arabic).

Other languages available on www.gimaitaly.com (PL, CZ, SK, HU, SI, HR, RS, RO, BG, RU).

Measurement range: 20-280 mmHg (pressure), 40-200 beat/minute (pulse)

Size: 138x94.5x62.5 mm

Weight: 402 g including batteries