

PREMIUM INFANT CUFF + 2 TUBE TPU BLADDER



|               |                |
|---------------|----------------|
| Product Code  | 32887          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|            |                |
|------------|----------------|
| UK-REP     | Yes            |
| CH-REP     | Yes            |
| RDM (NSIS) | 2512518        |
| CND        | C99            |
| EAN/UPC    | 8023279328875  |
| PARAF      | 987471404      |
| GTIN       | 08023279328875 |
| GMDN       | 34978          |

Description

PREMIUM CUFFS dark grey - nylon - TPU bladder - velcro fastening  
 Made of high quality nylon 210D.  
 Durability test up to 10,000 cycles.

Two tubes cuff for Minor, Yton line, Roma 2T, Dallas, Big Ben, Empire, Dayton...

- Model: Infant
- Size cm: 7.5x26
- Arm circumference: 13 - 19 cm