

ERKA SUPERB GREEN CUFF 2 TUBES - pediatric - for ERKA Vario



Product Code	32898
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
EAN/UPC	4250296150767
UMDNS	11-073

Description

ERKA GREEN CUFFS SUPERB RAPID 2 TUBES

- Child
- ERKA code: 025-26481
- Colour: green
- Size cm: 40x11.6
- Circumference cm: 20.5-28

Compatibility: Vario

Superb cuffs are easy to disinfect and clean, equipped with a waterproof surface. Cuffs are comfortable and consist of a cover and a removable inflation bladder. Super soft adhesive and fleece strap on the cuff. Made of durable patented German PU material. Optimal in blood pressure measurement because of pressure distribution and the highly precise values. PVC, DEHP and latex free. Available in different sizes with 1, 2 or twin tubes.

Made in Germany.