



Product Code	32901
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	No

RDM (NSIS)	2808935
CND	Z1203020501
EAN/UPC	8023279329018
UMDNS	16-157
PARAF	982171605

Description

ANDON B.P. MONITOR

- large easy to read LCD display, WHO blood pressure indicator
- 4 groups x 30 times with time and date stamp
- average of all saved measured blood pressure values for the last 7 days
- arrhythmia detection (HB)
- automatic switch-off

Multilingual box and manual: GB, FR, IT, ES, PT, DE

Cuff circumference: 22-30 cm

Technical Specifications

Measurement system: Oscillometric method

Measurement range:

Pressure: 0~300 mmHg \pm 3 mmHg

Pulse: 40~180 pulse/min \pm 5%

Power Source: 4 x 1.5 V AAA batteries - not included

Dimensions - weight: 107 x 80 x h 52 mm - 170 g (ex. batteries)

LCD size: 66 x 58.5 mm