



Product Code	32916
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No
RDM (NSIS)	2061057
CND	Z1203020501
EAN/UPC	8023279329162
PARAF	984594820

Description

GIMA BLUETOOTH BLOOD PRESSURE MONITOR

Bluetooth BPM with innovative tubeless technology and super slim design. It measures blood pressure and heart rate, detects irregular heartbeat and excessive body motion.

- Large LCD display 86x24 mm, blue mode display with backlight
- Comfortable cuff for arm circumference 22-42 cm (nylon outside + polyester inside in contact with skin)
- Date, time and low battery indication
- 2 groups of memory (2 users)

Multilingual box and manual: GB, FR, IT, ES.

Size: 131 x 73 x 29.4 mm

Technical Specifications

Measuring method: Oscillometric - continuous mode operation
 Measurement range: Rated cuff pressure 0~299 mmHg
 Measurement pressure: SYS: 60~230 mmHg, DIA: 40~130 mmHg / Pulse: 40-199 bpm
 Mechanical buttons: User 1 - User 2 - Query - Power
 Power supply: 3.7V - 1000 mHA built-in rechargeable li-polymer battery - battery life 2 years
 Charge power: up to 300 times
 AC Adaptor (included): Input 100-240V, 50-60Hz

COMPATIBILITY

For iOS devices: The operating system must be iOS 8 or higher, such as iPhone 4S, iPhone 5/5C/5S, iPhone 6/6 Plus and later models.
 For Android devices: The operating system must be 4.3 or higher.