

SMART GIMA AUTOMATIC WRIST BLOOD PRESSURE MONITOR



Product Code	32918
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	No

RDM (NSIS)	2868670
CND	Z1203020501
EAN/UPC	8023279329186
UMDNS	16-157
PARAF	927586368

Description

SMART AUTOMATIC WRIST B.P. MONITOR

Compact, fully logic, B.P.M. with innovative features:

- smart inflation technology for automatic pumping to proper pressure
- WHO blood pressure indicator
- IHB Irregular Heartbeat Detection
- average value of latest 3 readings
- 2x60 memories, date and time
- displays systolic, diastolic and pulse rate separately

Multilingual manual and box: GB, FR, IT, DE,ES, PT, PL, GR, Arabic

Display 31.9x39.6 mm LCD

Technical Specifications

Measuring method: oscillometric, automatic air inflation and measurement
Cuff pressure: 0-295 mmHg \pm 3 mmHg
Pulse rate range: 40~180 beats/min \pm 5% of reading
Range of measurement: 30-280 mmHg
Power supply: 2x"AAA" batteries 1.5 V
Dimensions - Weight: 84x64x28 mm - 110 g (excluding batteries)
Wrist cuff: applicable for wrist circumferences 14-19.5 cm