

SMART GIMA AUTOMATIC BLOOD PRESSURE MONITOR



Product Code	32921
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	No

RDM (NSIS)	2863334
CND	Z1203020501
EAN/UPC	8023279329216
UMDNS	16-157
PARAF	917437461

Description

SMART AUTOMATIC B.P. MONITOR with adult cuff

Compact, fully logic, multiuser, B.P.M. with innovative features:

- smart inflation technology for automatic pumping to proper pressure
- WHO blood pressure indicator
- IHB Irregular Heartbeat Detection
- average value of latest 3 readings
- 2x60 memories, date and time
- displays systolic, diastolic and pulse rate separately
- multilingual manual and box: GB, FR, IT, DE, ES, PT, PL, GR, Arabic

Display 62x43 mm LCD

Technical Specifications

Measuring method: oscillometric, automatic inflation and measurement

Measuring range:

- cuff pressure: 0-300 mmHg ± 3 mmHg
- pulse rate: 40~180 beats/min $\pm 5\%$ of reading
- systolic: 60-260 mmHg
- diastolic: 40-199 mmHg

Power supply: 4x"AA" batteries 1.5 V

Dimensions - Weight: 138x98x48 mm - 211 g (excluding batteries)

Arm cuff: for arm circumferences 22-30 cm, 30-42 cm or 42-48 cm