



Product Code	32926
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	Yes

RDM (NSIS)	409946
CND	Z1203020501
EAN/UPC	8023279329261
PARAF	926977416
GTIN	08023279329261
GMDN	45617

Description

Accurate and precise automatic digital blood pressure meters for smooth and fast measurements with innovative functions:

- IHB: irregular heartbeat detector
- WHO risk category indicator (classifies your measurements in 6 stages of risk, from optimal to severe hypertension)
- Proven clinically to meet European Society of Hypertension ESH protocol
- Easy to read, wide angle display
- 120 memories for 3 users: 3 zones of 40 memories each
- Start/stop unique button, with blue backlight
- Memory average: automatically average last 3 memories

Size: 74.5 x 74 x 30 mm

Weight: 95g

LCD size: 34.5 x 42 mm

Cuff size for arm circumferences: 13.5-19.5 cm

Technical Specifications

- Measuring method: oscillometric
- Measuring range:
 - Pressure: 0-300 mmHg
 - Pulse: 40-199 bpm
- Accuracy:
 - Pressure: ± 3 mmHg
 - Pulse: $\pm 5\%$
- Power supply: 2 x AAA alkaline
- Battery life: 250 measures