



Product Code	32936
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No
RDM (NSIS)	2556479
CND	Z1203020501
EAN/UPC	4015672101203

## Description

Professional B.P. monitor designed for clinical use in professional settings.

- Gold standard accuracy
- Validated according to international protocols:
  1. Association for the Advancement of Medical Instrumentation (AAMI)
  2. European Society of Hypertension
- International Protocol (ESH - IP)
- Intellisense technology
- Irregular pulse indicator
- 3 modes of operation: single mode, average mode, auscultation (manual mode)

Supplied with battery pack, adult cuff (22-32 cm), AC adapter.

Made in Japan.

## Technical Specifications

Method: Oscillometric

Measurement range: pressure: 0-299 mmHg  $\pm$  3 mmHg or 2%

pulse rate: 30-199 beats/min  $\pm$  5% of reading

Inflation: automatic inflation with pumping

Deflation: automatic deflation by electromagnetic control valve

Air release: automatic rapid air release by electromagnetic control valve

Pressure detection: electrostatic capacity semi-conductor pressure sensor

Power supply: AC adapter (230 VAC, 50 Hz, 20 VA)

Battery pack (4.8 VDC, 6W)

Electric shock protection method: Class II B type

Dimension: 139x203x131 mm (5.1/2"x8"x5.1/6")

Weight: approx. 910 g

Complies with: EN 1060, EN60601-1-1 - 2007