



Product Code	32964
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	I M

RDM (NSIS)	56383
CND	Z1302
EAN/UPC	8023279329643
PARAF	988245748

## Description

PHENOL FREE: completely free of phenols, it assures maximum safety to the user

Compatible with the following models of ECG:

- EDAN  
- SE-3
- DIMED  
- LTD405S, LTD405W, LTD415 Heartscr. 80G
- EDAN  
- SE-601
- INNOMED  
- HeartScreen IKO 80G-L, 80C-L

## Technical Specifications

Type: roll  
Width: 80 mm  
Length: 25 mm  
Sheets: -  
Marker: no  
Grid: orange