



Product Code	32965
Unit of sale	box of 20 pcs
Minimum order	1
Type	Medical device
Class	I M

RDM (NSIS)	56383
CND	Z1302
EAN/UPC	8023279329650
PARAF	988245635

Description

PHENOL FREE: completely free of phenols, it assures maximum safety to the user

Compatible with the following models of ECG:

- COMEN
- CM600
- DIMED
- LTD440
- LTD445/S
- EDAN
- SE-601
- FUKUDA ME - FUKUDA DENSHI
- OP222TE
- FX7202

Technical Specifications

Type: pack
Width: 110 mm
Length: 140 mm
Sheets: 143
Marker: yes
Grid: orange