



Product Code	32985
Unit of sale	box of 1
Minimum order	1
Type	Medical device
Class	I M

RDM (NSIS)	56383
CND	Z1302
EAN/UPC	8023279329858

Description

PHENOL FREE: completely free of phenols, it assures maximum safety to the user

Compatible with the following models of ECG:

ECG manufacturer: ESAOTE BIOMEDICA
Model: Archimed 4240, Archiwin

Technical Specifications

Type: pack
Width: 210 mm
Length: 300 mm
Sheets: 200
Marker: yes
Grid: orange