



Product Code	32988
Unit of sale	box of 20 pcs
Minimum order	1
Type	Medical device
Class	I M

RDM (NSIS)	56383
CND	Z1302
EAN/UPC	8023279329889

Description

PHENOL FREE: completely free of phenols, it assures maximum safety to the user

Compatible with the following models of ECG:

ECG manufacturer: FUKUDA ME - FUKUDA DENSHI
Model: 501A, FX121, FCP-500, OP18TE

ECG manufacturer: HEWLETT PACKARD
Model: 40453A, 43120

ECG manufacturer: KENZ-SUZUKEN
Model: 103, R50-30R, TEC 8300

Technical Specifications

Type: roll
Width: 50 mm
Length: 30 m
Sheets: -
Marker: no
Grid: orange