



Product Code	32996
Unit of sale	box of 20 pcs
Minimum order	1
Type	Medical device
Class	I M

RDM (NSIS)	56383
CND	Z1302
EAN/UPC	8023279329964

Description

PHENOL FREE: completely free of phenols, it assures maximum safety to the user

Compatible with the following models of ECG:

- ASPEL
 - Ascard B5 eco
 - Mr. Green
 - mint
- BTL
 - 08/MT
 - 08/MD
 - 08/MD3
- DIMED
 - LTD410
 - Heartscreen 60G
- GIMA
 - Cardiogima 1M
 - Cardiogima 3M
- INNOMED
 - HearthMirror IKO 3
 - HeartScreen 60 IKO
 - HeartScreen 60 G
- PROGETTI
 - EPG Basik-pocket

Technical Specifications

Type: roll
Width: 58 mm
Length: 25 m
Sheets: -
Marker: no
Grid: orange