



|               |                |
|---------------|----------------|
| Product Code  | 33019          |
| Unit of sale  | box of 1       |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I M            |

|            |               |
|------------|---------------|
| RDM (NSIS) | 56383         |
| CND        | Z1302         |
| EAN/UPC    | 8023279330199 |

## Description

PHENOL FREE: completely free of phenols, it assures maximum safety to the user

Compatible with the following models of ECG:

ECG manufacturer: MORTARA

Model: Eli 210, Eli 250, Eli 350, Eli 380

## Technical Specifications

Type: pack  
Width: 210 mm  
Length: 300 mm  
Sheets: 250  
Marker: yes  
Grid: orange