



| | |
|---------------|----------------|
| Product Code | 33028 |
| Unit of sale | box of 1 |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|------------|---------------|
| RDM (NSIS) | 2322675 |
| CND | Z1302 |
| EAN/UPC | 8023279330281 |

Description

PHENOL FREE: completely free of phenols, it assures maximum safety to the user

Compatible with the following models of ECG:

ECG manufacturer: NIHON KOHDEN

Model: 1500, 1550, 9110, 9130, 9132, 9320, 2250

Technical Specifications

Type: pack
Width: 210 mm
Length: 295 mm
Sheets: 170
Marker: yes
Grid: white