

# MEDIONE TRQ-2024 AUTOMATIC TOURNIQUET - 4 channels



Product Code	33110
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2789673
CND	Z12139006
EAN/UPC	8682754922245
GMDN	14074

## Description

### TRQ-2024 AUTOMATIC TOURNIQUET - 4 CHANNELS

Tourniquets used to temporarily stop or slow the blood flow in the patient's upper and lower extremities during surgical procedures. In addition, in the application of Regional Intravenous Anesthesia (IVRA), it prevents blood flow to that area after the blood of the extremity is drained. Designed for continuous use in adults and children. Adjustable pressure range is from 20 to 650 mmHg. Tourniquets have 4 connection sockets to use two cuffs. IVRA can be applied with a double cuff. Medical Device Class IIa.

#### Features:

- Automatic leakage test and system test
- Possibility to use 2 double or 4 single cuffs at the same time
- Safety, visual & audio alarm system
- 7-inch LCD touch screen
- Microprocessor Controlled
- Adjustable application time, pressure range and timing alarm
- Built-in rechargeable lithium-ion battery capacity 20 hours
- 4 tourniquet cuffs included
- Carrying handle

#### INTENDED USE OF TOURNIQUET SYSTEM

- Leg - arm fractures
- Mallet Finger therapy
- Arthroscopic Surgery
- Hammer toe
- Carpal Tunnel Syndrome
- Amputations
- Leg and arm implant insertion and removal
- Varicose vein
- Tendon correction
- Tumor excisions
- Finger, wrist, knee, elbow joint prosthesis
- Cyst excisions
- Implant insertion and removal for upper/lower extremities

## Technical Specifications

Cuff pressure range: 20 - 650 mmHg  $\pm$  5 mmHg  
Time range: 2-240 minutes  $\pm$  1 minute  
Nominal compressor value: 12 V, mask. 1.5 bar  
7" Touch Screen LCD 800x480 pixel  
Operating voltage: 100-240 V / 50-60 Hz 0.3A Max  
Rechargeable battery: Lithium battery, 14.4 V, 5.8Ah activates in case of power failure  
Size - Weight: 235x268xh 275 mm; 3.65 kg

## Standard accessories

1 reusable single cuff 12" - 30.5x7.6 cm - baby leg/child arm  
1 reusable single cuff 18" - 45.7x10.2 cm - child leg/adult arm  
1 reusable single cuff 30" - 76.2x12.7 cm - adult  
1 reusable IVRA double cuff 20" - 50.8x15.2 cm - adult  
2 blue patient connecting tubes - 3 m  
2 red patient connecting tubes - 4 m  
User manual GB, FR, IT, ES. TR available on [gimaitaly.com](http://gimaitaly.com)  
Power cable