

# NO PINCH SINGLE CONE CUFF 86 x 11 cm - thigh L - red - autoclavable



Product Code	33137
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

GMDN	17230
------	-------

## Description

CUFFS FOR DTS-2000 AND DTS-3000 PNEUMATIC TOURNIQUETS

### SINGLE CUFF

- Autoclavable: Silicone rubber bladder
- 86x11 cm - red
- Size Arm: -
- Size Thigh: L

- no pinching, no damage to the skin
- lower pressure compared to normal cuff (50 mmHg lower)
- no leaking (main reason for lower pressure)