

# NO PINCH SINGLE CUFF 61 x 9 cm - arm L or thigh S - red - spare



Product Code	33141
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
EAN/UPC	8023279331417
GMDN	17230

## Description

CUFFS FOR DTS-2000 AND DTS-3000 PNEUMATIC TOURNIQUETS

### SINGLE CUFF

- Reusable: EtCO2 sterilization
- 61x9 cm - red
- Size Arm: L
- Size Thigh: S

- no pinching, no damage to the skin
- lower pressure compared to normal cuff (50 mmHg lower)
- no leaking (main reason for lower pressure)