

# NO PINCH DOUBLE CUFF 57 x 10 cm - arm M - blue - spare



|               |                |
|---------------|----------------|
| Product Code  | 33145          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|        |    |
|--------|----|
| UK-REP | No |
| CH-REP | No |

|         |               |
|---------|---------------|
| EAN/UPC | 8023279331455 |
| GMDN    | 17230         |

## Description

### CUFFS FOR DTS-2000 AND DTS-3000 PNEUMATIC TOURNIQUETS

#### DOUBLE CUFF

- Reusable: EtCO<sub>2</sub> sterilization
- 57x10 cm - blue
- Size Arm: M
- Size Thigh: -

- no pinching, no damage to the skin
- lower pressure compared to normal cuff (50 mmHg lower)
- no leaking (main reason for lower pressure)