

# NO PINCH DOUBLE CUFF 80 x 15 cm - thigh M - red - spare



Product Code	33146
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
EAN/UPC	8023279331462
GMDN	17230

## Description

### CUFFS FOR DTS-2000 AND DTS-3000 PNEUMATIC TOURNIQUETS

#### DOUBLE CUFF

- Reusable: EtCO2 sterilization
- 80x15 cm - red
- Size Arm: -
- Size Thigh: M

- no pinching, no damage to the skin
- lower pressure compared to normal cuff (50 mmHg lower)
- no leaking (main reason for lower pressure)