



Product Code	33150
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	1874053
CND	Z12139006
EAN/UPC	4045396007499

Description

Reliable pneumatic tourniquet used during amputations to restrict blood circulation in the limbs. It has a robust and very accurate manometer to provide exact measurements and a particularly long service life.

- Sturdy metal pump to produce pressure
- Air in the cuff can escape by opening the release valve
- Chromium-plated metal manometer with easy-to-read aluminium scale up to 700 mm Hg

Pneumatic tourniquet includes 2 cuffs (upper arm cuff 57 x 9 cm and thigh cuff 96 x 13 cm) with latex free bladder.