

SAM PRO TRAINER for Semi-Automatic Rescue Sam AED Defibrillator- Italian



Product Code	33171
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279331714

Description

Designed to effectively train emergency responders to the use of RESCUE SAM, thanks to simple, clear, spoken and visual instructions.

5 pre-programmed scenarios selected directly by the remote control.

Totally safe real shock simulation.

Supplied with remote control, training pads, battery charger.