



Product Code	33176
Unit of sale	box of 6 pcs
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No
EAN/UPC	8023279331769
PARAF	977813500

Description

Special dermatograph pencil to mark the skin for identification of areas to be treated.
Made in cedar wood with wax lead.

Colour: white (6 pieces)

Size: 175 x Ø 7 mm.

Suitable for general surgery, dermatology, orthopaedy, radiology...

Pencils are sold non sharpened in box of 6 with 1 pencil sharpener included.

Multilanguage box and manual: GB, FR, IT, DE, ES, PT, GR, PL, SE, RO, HU, Arab.