


GIMA DERMATOGRAPH PENCILS - mix colours



Product Code	33178
Unit of sale	box of 6 pcs
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No

EAN/UPC	8023279331783
PARAF	977813524

Description

Special dermatograph pencil to mark the skin for identification of areas to be treated.
Made in cedar wood with wax lead.

Mixed colors: red, blue, white (2 pieces each)
Size: 175 x Ø 7 mm.

Suitable for general surgery, dermatology, orthopaedy, radiology...

Pencils are sold non sharpened in box of 6 with 1 pencil sharpener included.

Multilanguage box and manual: GB, FR, IT, DE, ES, PT, GR, PL, SE, RO, HU, Arab.