



Product Code	33179
Unit of sale	box of 10 pcs
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No

EAN/UPC	8023279331790
PARAF	977813536

Description

SKIN MARKER - dual tips gentian violet - box of 10

Double tip skin marker with polypropylene cap and body.
Polyester fiber tips.

- Ultra fine tip utility marker to write on adhesive labels and infusion bags.
- Regular tip surgical skin marker for use on patient skin.

Washable with water and soap.
Size 138xØ 12.5 mm.

Multilingual box and manual: GB, FR, IT, DE, ES, PT, GR, PL, SE, RO, HU, Arabic.