

SURGICAL SKIN MARKER - single tip 1.0 mm - sterile



Product Code	33182
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	IS

RDM (NSIS)	2347171
CND	V9004
EAN/UPC	8023279331820

Description

SURGICAL SKIN MARKER - single tip - sterile

- Single tip marker ideal for plastic surgery.
- Hypoallergenic and violet non-toxic ink.
- Marking can be removed by alcohol and liquid soap.
- Intended for single use.
- Includes a separate six inch, flexible ruler calibrated in both inches and centimeters.
- Individually packed.

Technical Specifications

- Tip size: 1.0 mm
- Diameter: 9.5 mm
- Length: 167 mm
- Ink colour: violet
- Washable: no