

# SURGICAL SKIN MARKER - double tip 0.5 and 1.0 mm - sterile



Product Code	33183
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	I S

RDM (NSIS)	2347171
CND	V9004
EAN/UPC	8023279331837

## Description

SURGICAL SKIN MARKER - double tip - sterile - box of 100

Double tip marker ideal for plastic surgery. Hypoallergenic and violet non-toxic ink. Marking can be removed by alcohol and liquid soap. Intended for single use.

Includes a separate six inch, flexible ruler calibrated in both inches and centimeters. Individually packed.

Multilingual box and instructions: GB, FR, IT, ES, PT, DE, NL, DK, SE, NO, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, Arabic.

## Technical Specifications

Tip size: fine tip: 0.5mm / thick tip: 1.0mm  
Diameter: 10.2 mm  
Length: 138.5 mm  
Ink colour: violet  
Washable: no