

SURGICAL SKIN MARKER - single tip 1.0 mm - sterile



Product Code	33184
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	I S

UK-REP	No
CH-REP	No

RDM (NSIS)	2347171
CND	V9004
EAN/UPC	8023279331844
PARAF	965138656

Description

SURGICAL SKIN MARKER - single tip - sterile - box of 100

Regular tip.

Single tip marker ideal for plastic surgery. Hypoallergenic and violet non-toxic ink. Marking can be removed by alcohol and liquid soap. Intended for single use.

Includes a separate six inch, flexible ruler calibrated in both inches and centimeters. Individually packed.

Multilingual box and instructions: GB, FR, IT, ES, PT, DE, NL, DK, SE, NO, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, Arabic.

Technical Specifications

Tip size: 1.0 mm
Diameter: 9.5 mm
Length: 167 mm
Ink colour: violet
Washable: no