



Product Code	33206
Unit of sale	box of 500 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	60766
CND	C020501
PARAF	977813587
GTIN	05707480028047
GMDN	35035

Description

- AMBU BLUE SENSOR VL
- highly conductive wet gel
 - combination of instant and long-term adhesives
 - breathable microporous material
 - unique offset connector
 - high quality Ag/AgCl sensor

BlueSensor is the world’s most time and cost efficient range of electrodes that delivers unmatched readings, user convenience and patient comfort. BlueSensor ensures that you don’t get your workflow disrupted by false alarms or have to waste time resolving problems with electrode adhesion.

Exact monitoring clear picture Blue Sensor electrodes give a fast and reliable trace without the need for skin preparation, maximizing comfort for the patient. We know that fast, exact monitoring that ensures a clear picture of the patients’ heartbeat is essential to accurate diagnoses.

Technical Specifications

- AMBU code: VL-00-A/25
- Description: BLUE SENSOR VL
- Application: Long-term
- Gel type: Wet
- Backing material: Microporous
- Size(mm): 72 x 68
- Shape: Round with offset
- Shelf life* (month): 24
- Radiolucent: NO
- Connector type: 4 mm
- Minimum order: Box of 500 pcs.
- Packaging: 20 pouches of 25 pcs.

*Shelf life means maximum life from production date. Average Gima shelf life is 50 to 90% of maximum life.