



Product Code	33222
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2962141
CND	Z12050301
EAN/UPC	6945040100485
UMDNS	11407
PARAF	977813649

Description

- 7" TFT, colour LCD touch screen
- 12 leads simultaneous display
- Digital isolation technology and signal processing solution, digital filter
- Patient management, name, age and id, detail analysis report
- Built-in lithium rechargeable battery
- Thermal printing system provides more than 10 printing modes
- Can print 90 minutes and print out 150 ECG continually
- ECG-Sync software, USB connection to PC and built-in memory can store more than 1000 cases
- Auto-measurement, auto-analysis and auto-interpretation
- Internal software: GB, FR, IT, ES, DE, TR, PL, PT, SI, RU, KZ, CN.
- PC software ECG Synchronous: GB, IT, FR, DE, PL, SK, TR.
- Manual: GB, IT, FR.

USB connection allows real-time transmission to PC, remote data analysis and diagnosis, management and printing through the PC.

Technical Specifications

ECG instrument Input

- Safety Type: Class I, type CF applied part, which has defibrillation-proof applied part
- Lead Standard 12 leads
- Input way: Floating and defibrillation protection
- Patient leak current: <math><10\mu\text{A}</math>
- Sampling rate: 1000Hz
- A/D conversion: 12bit
- Frequency response 0.05Hz~150Hz (-3dB, + 0.4dB)
- Time constant: >3.2s
- CMRR: >60dB, >100dB(Add filter)
- Enduring polarization voltage: $\pm 500\text{mV}$
- Noise level: <math><15\mu\text{Vp-p}</math>

Dispose Wave Data

- EMG interference filter: 25Hz/35Hz (-3dB)
- AC filter 50Hz/60Hz (-20dB)

Waveform output

- Recording way: Thermal printing system (32spots/mm)
- Recording paper: 110mm(W) \times 20m(L)
- Sensitivity: 2.5, 5, 10, 20, 40 mm/mV
- Paper speed:
 - Auto-record: 25 mm/s, 50 mm/s
 - Rhythm record: 25 mm/s, 50 mm/s
 - Manual-record: 5, 6.25, 10, 12.5, 25, 50, 100 mm/s

Power supply

- AC: 110-220V, 50/60 Hz
- DC: Lithium battery 7.4V, 3500 mAh

Fuse

- 2 pcs Ø5×20mm AC time lag; T250mA/250V(Power Supply:220V)

Form

- Size 315 x 215 x 77 mm
- Net Weight 1.7 kg

Standard accessories

- ECG cable - 1 set
- Limb electrode clamp - 1 set
- chest electrodes - 1 set
- Ground wire
- Power connection
- Paper roll
- USB connection and PC software