



Product Code	33238
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	No

RDM (NSIS)	2573615
CND	Z1203020202
EAN/UPC	6934440102175
PARAF	950453922
GMDN	30004

Description

DuoEK™ S PORTABLE ECG MONITOR

Portable ECG monitor for instant ECG analysis and measurement.
 Compact and lightweight, it easily performs medical-grade ECG examinations anytime, anywhere.
 3 measurement leads: hand-hand, hand-leg and chest.
 - OLED screen for clear ECG display
 - one-touch measurement
 - 500 measurements on a single charge

Suitable for people with high stress, irregular schedules, heart disease and heart discomfort (palpitations, panic, dizziness, chest tightness, etc.).

Free ViHealth App on iOS and Android
 Automatically upload data and share PDF report with doctor.

Supplied with charger and multilingual manual (GB, FR , IT, ES, DE).

Technical Specifications

ECG length: 30 seconds to 5 minutes recording
 Heart rate range: 30 to 250 bpm
 Battery: rechargeable lithium polymer battery > 500 checks
 Charging time: 2-3 hours
 Display: 0.96" OLED display
 Connection: micro-USB connector
 Wireless: Bluetooth 4.0 BLE
 Built-in memory:
 - 10 sessions, if ECG length < 1 minute
 - 4 sessions, if ECG length is 5 minutes
 App storage: unlimited
 Size - weight: 92x32x8.2 mm - 26.8 g