



Product Code	33256
Unit of sale	box of 2 pcs
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No
EAN/UPC	8683229603089

## Description

### KONIX KAVITATION GEL

Cavitation gel - bag 5 l  
Minimum order: box of 2  
Includes 1 empty bottle 250 ml and 1 nozzle for easy refill

Cavitation gel reduces the pain that patients experience during depilation and makes the device last longer by cooling the device head and window.

Thanks to the vacuum technology, its active substances and its high viscosity preserve probe from damages.

It is not affected by the body salt during the use, ensures an easy application.

Technical properties: it is hypoallergenic, does not irritate skin, is water soluble and easily cleanable.

Oily substances, formaldehyde and salt free.

No toxic effects and odourless.