



Product Code	33261
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

CH-REP	Yes
--------	-----

RDM (NSIS)	1980609
CND	Z12050301
EAN/UPC	8023279332612
PARAF	977813740
GMDN	62641

## Description

Designed for general practitioners or for domestic users who suffer from cardiovascular diseases such as hypertension, coronary disorders, giving the possibility of a daily ECG test. Records can be saved and used for clinical applications.

4 measurement modes: palm, leg, chest or lead wire traditional measurement.

Main features:

- small, portable and easy to use
- 30 second ECG waveform and measurement results
- 1200 ECG waveform can be saved
- 10 hours ECG continuous measurement data storage
- analysis and report of 17 types of ECG waveform
- USB 2.0 interface with optional data cable and ECG software
- traditional measurement with optional ECG cable

Supplied with multilanguage manual and box (GB, FR, IT, ES).

## Technical Specifications

- ECG bandwidth: 0.5 Hz - 40 Hz
- Heart rate measuring range: 30 - 240 bpm  $\pm$  2 bpm
- LCD back light display: 240 x 160 dot matrix
- Display size: 57.6 x 40mm
- Power: 2 x AAA batteries
- Size 125 x 70 x 21.5 mm
- Weight: 106 g