

EEG PASTE - 400 g



Product Code	33262
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2918897
CND	V9015
EAN/UPC	8683229601153
GMDN	11425

Description

EEG PASTE - 400 g

Used with EEG electrodes in EEG application for brain mapping. Thanks to the active substances in its content, it allows the electrode to cling to the strand of hair easily and provides conductivity.

It is hypoallergenic, does not irritate skin, soluble in water and can be cleaned easily.

Does not contain oil and oily substances.

It has no toxic effects and odour.