



Product Code	33267
Unit of sale	box of 2 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2918897
CND	V9015
EAN/UPC	8683229603072
GMDN	11425

Description

ECG GEL - 5 l - box of 2 bags
(Includes 1 empty bottle 250 ml and 1 nozzle for easy refill)

Suitable for non pre-gelled electrodes, ECG latches and pendant switches in any ECG and exertional applications. Allows ECG waves clear transfer to equipment. It is hypoallergenic, does not irritate skin, does not contain oil and oily substances, does not contain formaldehyde and salt, has no toxic effects and is odourless.