



Product Code	33269
Unit of sale	box of 48 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No
RDM (NSIS)	2958448
CND	V9099
EAN/UPC	8681715349695
PARAF	977813753
GMDN	58735

## Description

### GIMA ULTRASOUND GEL

Sterile ultrasound gel - sachet - 20ml - transparent  
Minimum order: box of 48

Hypoallergenic, water soluble gel suitable for ultrasonography, doppler and defibrillator applications. It ensures the ultrasound waves come to the equipment display in a clear and uninterrupted way. Salt, formaldehyde, fragrance free.