

ULTRASOUND GEL - tube 60 ml - transparent



Product Code	33270
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2971718
CND	V9099
EAN/UPC	8023279332704
PARAF	977813765
GMDN	15321

Description

GIMA ULTRASOUND GEL

Ultrasound gel - tube 60 ml - transparent
Minimum order: box of 12

Hypoallergenic, water soluble gel suitable for ultrasonography, doppler and defibrillator applications. It ensures the ultrasound waves come to the equipment display in a clear and uninterrupted way. Salt, formaldehyde, fragrance free.