

ULTRASOUND GEL - tank 5 l - blue



Product Code	33276
Unit of sale	box of 2 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2450834
CND	V9099
EAN/UPC	8023279332766
PARAF	977813815
GTIN	08023279332766
GMDN	15321

Description

GIMA ULTRASOUND GEL

Ultrasound gel - tank 5 l - blue
Minimum order: box of 2
Includes 2 empty bottles 250 ml

Hypoallergenic, water soluble gel suitable for ultrasonography, doppler and defibrillator applications. It ensures the ultrasound waves come to the equipment display in a clear and uninterrupted way. Salt, formaldehyde, fragrance free.