

# ULTRASOUND GEL - tank 5 l - transparent



Product Code	33277
Unit of sale	box of 2 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	2450834
CND	V9099
EAN/UPC	8023279332773
PARAF	977813827
GTIN	08023279332773
GMDN	15321

## Description

### GIMA ULTRASOUND GEL

Ultrasound gel - tank 5 l - transparent  
Minimum order: box of 2  
Includes 2 empty bottles 250 ml

Hypoallergenic, water soluble gel suitable for ultrasonography, doppler and defibrillator applications. It ensures the ultrasound waves come to the equipment display in a clear and uninterrupted way. Salt, formaldehyde, fragrance free.