

# ULTRASOUND GEL - bag 5 l - blue



Product Code	33286
Unit of sale	box of 2 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2450834
CND	V9099
EAN/UPC	8023279332865
GTIN	08023279332865
GMDN	15321

## Description

### GIMA ULTRASOUND GEL

Ultrasound gel - bag 5 l - blue

Minimum order: box of 2

Includes 1 empty bottle 250 ml and 1 nozzle

Hypoallergenic, water soluble gel suitable for ultrasonography, doppler and defibrillator applications. It ensures the ultrasound waves come to the equipment display in a clear and uninterrupted way. Salt, formaldehyde, fragrance free.