

ULTRASOUND GEL - bottle 1 l - blue



Product Code	33288
Unit of sale	box of 18 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2450834
CND	V9099
EAN/UPC	8023279332889
GTIN	08023279332889
GMDN	15321

Description

GIMA ULTRASOUND GEL

Ultrasound gel - bottle 1 l - blue
Minimum order: box of 18

Hypoallergenic, water soluble gel suitable for ultrasonography, doppler and defibrillator applications. It ensures the ultrasound waves come to the equipment display in a clear and uninterrupted way. Salt, formaldehyde, fragrance free.