

3M RED DOT 2228 FOAM ELECTRODES - 3.71x3.25 cm
- 20 bags of 50 pcs



Product Code	33290
Unit of sale	box of 1000 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2238458
CND	C020501
GTIN	50707387789522

Description

3M™ RED DOT™ ELECTRODES

Foam monitoring electrodes

- 3M Code: 2228
- 3M™ Red Dot™ Electrodes - 3.71 x 3.25 cm
- Box of 1,000 pieces (20 bags of 50 pieces)

3M 2228 adhesive foam-tape backing and proprietary sticky gel ensure excellent skin contact. Radiolucent. 3-days wearing time.