

# 3M RED DOT 2228BA FOAM ELECTRODES - 3.71x3.25 cm 20 bags of 40 pcs



Product Code	33291
Unit of sale	box of 800 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2238461
CND	C020501
GTIN	54046719263511

## Description

### 3M™ RED DOT™ ELECTRODES

Foam monitoring electrodes

- 3M Code: 2228BA
- 3M™ Red Dot™ Electrodes with 4 mm adaptor - 3.71 x 3.25 cm
- Box of 800 pieces (20 bags of 40 pieces)

3M 2228BA adhesive foam-tape backing and proprietary sticky gel ensure excellent skin contact. Thickness allows quick application build-up and padding over irregular connections.