

3M RED DOT 2244 FOAM ELECTRODES - 4x3.25 cm



Product Code	33292
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2238475
CND	C020501
GTIN	10707387789449

Description

3M™ RED DOT™ ELECTRODES

Foam monitoring electrodes

- 3M Code: 2244
- 3M™ Red Dot™ Electrodes - 4 x 3.25 cm
- Bag of 50 pieces

3M 2244 adhesive foam-tape backing and proprietary sticky gel ensure excellent skin contact. Radiolucent. 3-days wearing time.