

# 3M RED DOT 2570 FOAM ELECTRODES - 4x3.5 cm



Product Code	33293
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2240048
CND	C020501
GTIN	10707387789319

## Description

### 3M™ RED DOT™ ELECTRODES

Foam monitoring electrodes

- 3M Code: 2570
- Size: 4 x 3.5 cm
- Electrode with foam tape and sticky gel - radiolucent. 5-days wearing time.