

# 3M RED DOT 2238 NON WOVEN ELECTRODES soft cloth - 6 cm diam - 20 bags of 50 pcs



Product Code	33296
Unit of sale	box of 1000 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2120052
CND	C020501
GTIN	50707387788846

## Description

### 3M™ RED DOT™ ELECTRODES

Non woven monitoring electrodes

- 3M Code: 2238
- 3M™ Red Dot™ - Ø 6 cm - soft cloth
- Box of 1000 pieces (20 bags of 50 pieces)

3M 2238 soft cloth monitoring electrode with solid gel and lift tab. Breathable, flexible and conformable, gentle to the skin for patient comfort. Radiolucent. 3-days wearing.