

# 3M RED DOT 2239 NON WOVEN ELECTRODES

micropore - 6 cm diam - 20 bags of 50 pcs



Product Code	33297
Unit of sale	box of 1000 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2238469
CND	C020501
PARAF	987474590
GTIN	50707387788853

## Description

### 3M™ RED DOT™ ELECTRODES

Non woven monitoring electrodes

- 3M Code: 2239
- Size: Ø 6 cm
- Monitoring electrode with micropore tape backing. 3-days wearing time.
- 20 box of 50 pcs.