

# 3M RED DOT 2248-50 ELECTRODES - 4.5 cm diam



Product Code	33299
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2238477
CND	C020501
PARAF	989251069
GTIN	10707387789258

## Description

3M™ RED DOT™ ELECTRODES Pediatric

Non woven monitoring electrodes

- 3M Code: 2248-50
- Size: Ø 4.5 cm
- Pediatric monitoring electrode with micropore tape backing.